

Product Instruction

Product name:AI Fever Screening SystemModel:AWS-60-E12

Appearance:




AI Detector(K03)Intelligent Terminal(B39-YT)

Wiring operation:



Power-over-EthernetAdapter12VPower Adapter



- ① Connect Power-over-Ethernet Adapter
- ② Connect 12V Power Adapter
- ③ Connect AI Detector's RJ45 network port
- ④ Connect AI Detector's USB port
- ⑤ Connect mouse
- ⑥ Connect the monitor with an HDMI cable
- ⑦ Connect 12V Power Adapter

AI Fever Screening System hardware component (AWS-60-E12):

AI Detector (K03) + Intelligent Terminal (B39-YT) + Monitor (Optional)

Hardware specification:

			
Language	Multi-language		
Processor	RK3399, six-core, dual-core A72 + quad-core A53, @ 1.8GHz		
RAM	DDR3 4GB		
ROM	EMMC 16GB		Maximum support 128GB
HDMI output	1, support 4K/1080P		
camera	N/A		
Memory card slot	1, TF Card		Maximum support 128GB
Headphone slot	1, AUDIO		
Microphone	N/A		
USB	1USB 2.0、1USB3.0、1Type-C		
Network	2.4G WiFi		
	10/100M/1000M Ethernet		Support dual network ports
Characteristics			
Video format	Support mainstream video formats such as wmv, avi, flv, rm, rmvb, mpeg, ts, mp4		
Video decoding	Support 1080P, 4K, H.265 hard decoding		
Playback mode	Support multiple playback modes such as loop, timing, insert		
Audio format	Support audio formats such as MP3		
Image Format	Support BMP, JPEG, PNG, GIF, etc.		
Power Adapter	DC12V, 2A		
RTC	Support real-time clock		
Scheduled power on/off	Support		
Boot method	Connect the power to boot		
Software support			
Operating system	Android 8.1		
Startup screen	Android default		
Installation	APK Installer		

System Management	System settings		
Time	Global time		
System Upgrade	Support TF card / USB		
Input	Standard Android keyboard, optional third-party input methods supported (Chinese, Korean, Japanese, etc.)		
App	Standard Android Default Apps		
Exterior			
Color	Black, metal		



Temperature measurement	
Measurement range	32°C~42°C
Accuracy	± 0.5 °C (target temperature: 32 °C ~ 42 °C, fluctuate ± 0.3 °C under certain conditions)
Measure time	Less than 1 second
Measure distance	Within 2 meters, the optimal distance is 1 ~ 1.5 meters
Interface	
Thermometer	USB + RJ45 (with power supply)
Facial recognition camera	USB
External interfaces of terminal	RJ45、WIFI、HDMI
Environmental adaptability	
Ambient temperature	The best ambient temperature is 16 ~ 32 °C
Storage temperature	-20~60°C
Working humidity	< 90% (non-condensing)
Face recognition	Local recognition
Face library size	10000 faces database as standard

Product description:


This is a complete solution for face recognition + temperature screening + big data statistics. The AI detector can quickly capture the face information and temperature information. Then the intelligent terminal calculates the body temperature, identifies the captured face and displays the result. When the body temperature exceeds the normal body temperature set, an alarm will be issued to remind the relevant personnel and medical staffs, and the abnormal temperature data will be uploaded to the server for data statistics at the same time. In conclusion, users can check all members' temperatures and prevent the disease outbreak through the product.

UI description



5 Modules

① Statistics: View temperature statistics of the day or days.



② Personnel: Register, query, edit and delete personnel

1. Enter the name, query the personnel's temperature.

2. Click "≡" to add, batch import or batch delete personnel.

3. Note: The server version registers personnel through the server.







Photo Registration Batch Import



- **Photo registration:** Take a photo or choose an image to register.
- **Batch Import:** Upload a folder with member photos named as "Name_Number".
- **Batch Deletion:** Select the personnels to be deleted and click Delete.

③Record: Select time to view historical temperatures or to export that. Clear data.



④Settings

- **Recognition mode:** Switch single mode and multiple mode.
- **Device information:** View the device information and check for update.
- **Online authorization:** Default "authorized".
- **Advance:** If the temperature deviation is too large due to environmental influence, please calibrate under the guide of a technician.

⑤Log: If the system fails, please upload the logs, and the technician will analyze the cause to solve the problem

FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

The distance between user and device should be no less than 20cm.