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
Manuel d'instructions

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
PULSE OXIMETER

Instruction manual

patients' SpO₂ and pulse rate, please clip the thick finger such as thumb and middle finger deeply enough into the probe.

(29)The update period of data is less than 5 seconds, when the measured data is beyond the highest or lowest limit, the device would start alarming automatically on the premise of the alarming function is on.

(34)The device may not work for all patients. If you are unable to achieve stable readings,discontinue use.

(35)A flexible circuit connects the two parts of the device. Do not twist or pull on the connection.

(36)Take out batteries if the equipment is not used for a long time.

(37)Properly carry the equipment and prevent it from falling, collision, strong vibration or other mechanical external force damage.

(14)As the measur is taken on the basis of arteriole pulse,substantial pulsating blood flow of subject is required. For a subject with weak pulse due to shock, low ambient/body temperature, major bleeding, or use of vascular contracting drug, the SpO₂ waveform (PLETH) will decrease. In this case, the measurement will be more sensitive to interference.

4.BLUETOOTH AND APP

4.1 Download instructions

1.Scan the code and install the app, as shown in the figure below.(If you scan the QR code and click install, you cannot download it, You can choose to open it in the browser and download it.)

4.2 Operation interface description

1.Click the "connect to device" button.Then switch to the measurement interface.

6.Get the information directly from screen display.

7.The button has two functions.When the device is in standby mode, pressing the button can exit it;When the device is in operation status, pressing the button long can change brightness of the screen.

8.The device could change display direction according to the handing direction.

8.REPAIRING AND MAINTENANCE

Please change the batteries when the low-voltage displayed on the screen.

Please clean the surface of the device before using,Wipe the device with medical alcohol first and then let it dry in air or clean it by dry clean fabric.

Using the alcohol to disinfect the product after use, prevent from cross infection for next time use.

Please take out the batteries if the oximeter is not in use for a long time.

The best storage environment of the device is40°C to 60°C ambient temperature and not higher than 95% relative humidity.

Users are advised to calibrate the device termly (oraccording to the calibrating program of hospital). It also can be performed at the state-appointed agent or just contact us for calibration.

INSTRUCTIONS TO USER

Dear Users, thank you very much for purchasing our product.

The Manual is written for the current Pulse Oximeter. in case of modifications and software upgrades, the information contained in this document is subject to change without notice.

The Manual describes. in accordance with the Pulse Oximeter's features and requirements, main structure, functions, specifications, correct methods for transportation, installation, usage, operation, repair, maintenance and storage, etc.as well as the safety procedures to protect both the user and equipment. Refer to the respective chapters for details.

Please read the Manual very carefully before using this equipment. These instructions describe the operating procedures to follow these instructions can cause measuring abnormally, equipment damage and personal injury.

The manufacturer is not responsible for the safety, reliability and performance issues and any monitoring abnormality, personal injury and equipment damage due to users negligence of the operation instructions. The manufacturer's warranty service does not cover such faults.

Owing to the forthcoming renovation, the specific products you received may not be totally in accordance with the description of this User Manual. We would sincerely regret for that.

2.OVERVIEW

The pulse oxygen saturation is the percentage of HbO₂ in the total Hb in the blood,so-called the O₂ concentration in the blood. It is an important bio-parameter for the respiration. For the purpose of measuring the SpO₂ more easily and accurately, our company developed the Pulse Oximeter. At the same time, the device can measure the pulse rate simultaneously.

The Pulse Oximeter features in small volume, low power consumption, convenient operation and being portable. It is only necessary for tester to put one of his fingers into a fingertip photoelectric sensor for diagnosis, and a display screen will directly show measured value of Hemoglobin Saturation.

2.1 Features


Operation of the product is simple and convenient.

The product is small in volume, light in weighttotal weight is about 28g including batteries land convenient in carrying.



Power consumption of the product is low and the two originally equipped AAA batteries can be operated continuously for 20 hours.

The product will enter standby mode when no signal is in the product within 5 seconds.

Display direction can be changed, easy to view.

2.Click the "History" button.Then switch to the historical data view interface.

9.TROUBLESHOOTING

Trouble	Possible Reason	
The SpO ₂ and Pulse Rate are not displayed normally.	1.The finger is not properly positioned. 2.The tested SpO ₂ is too low to be detected.	Place the finger properly and try again.
The SpO ₂ and Pulse Rate are not displayed stably.	1.The finger is not placed inside deep enough. 2.The finger is shaking or the tester is moving.	1.Place the finger properly and try again. 2. Let the tester keep calm.
The device cannot be turned on.	1.The batteries are drained or almost drained. 2. The batteries are not inserted properly. 3.The malfunction of the device.	1. Change batteries. 2. Reinsert batteries. 3. Please contact the local service center.
The display is off suddenly.	1.The product will enter standby mode when no signal is in the product within 5 seconds. 2.The batteries are almost drained.	1. Normal. 2. Change batteries.

1.SAFETY INFORMATION

1.1 Warning

- (1)Only qualified personnel designated by the manufacturer is allowed to maintain the equipment. Users themselves are not allowed to repair the equipment.
- (2)Explosion hazard: not use the equipment in environment with inflammable substances such as anesthetics.
- (3)Not use the equipment in MRI and CT examination.
- (4)When used with electro-surgical equipment together, the user shall ensure the safety of the monitored patient.
- (5)Users who are allergic to silicone rubber is forbidden to use this equipment.
- (6)The equipment, accessories and packaging (batteries, plastic bags, foam and cartons, etc.) shall be scrapped in compliance with local laws and regulations.
- (7)Forbidden to use the equipment once the equipment is found damaged, or have material deterioration.
- (8)Do not use the equipment under the application conditions beyond its declared specification scope.
- (9)When using the equipment for long-term continuous monitoring, check and replace measured parts at least every 2 hours. It is necessary for more frequent examinations in some scenarios, such as perfusion disorders and skin sensitivity. Prolonged monitoring may increase unpredictable skin changes, such as allergies,

2.2 Major Applications and Scope of Application

The Pulse Oximeter can be used to measure human Hemoglobin Saturation and pulse rate through finger, and indicate the pulse intensity by the bar-display. The product is suitable for use in family, Oxygen Bar, social medical organizations and also the measure of saturation oxygen and pulse rate.

2.3 Environment Requirements

- Storage Environment
- a)Temperature: -40°C~+60°C
 - b)Relative humidity: <95%
 - c)Atmospheric pressure: 500 hPa~1060 hPa
- Operating Environment
- a)Temperature: 10°C~40°C
 - b)Relative humidity: <75%
 - c)Atmospheric pressure: 700 hPa~1060 hPa

3.PRINCIPLE AND CAUTION

3.1 Principle of Measurement

Principle of the Oximeter is as follows: An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive Hemoglobin (Hb) and Oxylhemoglobin (HbO₂) in glow & near-infrared zones. Operation principle of the instrument is: Photoelectric

reddening, blistering, or compressive necrosis.

(10)For the special patients, there should be a more prudent inspecting in the placing process.The device cannot be clipped on the edema and tender tissue.

(11)The light (the infrared is invisible) emitted from the device is harmful to the eyes, so theuser and the maintenance man should not stare at the light.

(12)Testee cannot use enamel or other makeup.

(13)Testee's fingernail cannot be too long.

(14)Ensure that the environment in which the device is operated is not subject to any sources of strong electromagnetic interference, such as radio transmitters, mobile telephones, etc. Keep them far away High level electromagnetic radiation emitted from such devices may greatly affect the monitor performance.


(15)Be careful with the use of the lanyard cord, Improper use of the lanyard cord will cause device damage not covered under the manufacturer's warranty. Swinging the device by the lanyard cord will void the warranty. Please do not use lanyard cord if allergic to lanyard cord.

(16)Please choose the accessories and probe which are approved or manufactured by the manufacturer, or else it may damage the device.

(17)Please don't measure this device with functional tester for the device's related information.

(18)When uploading, do keep the patient form touching the USB port. The computer used when uploading must be in accordance with EN60950.We recommended to use Lenovo Qitian M320,etc.in addition, when the date line

Oxylhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning Recording Technology, so that two beams of different wavelength lights can be focused onto human nail tip through perspective clamp fingertype-sensor. Then measured signal can be obtained by a photo-sensitive element,information acquired through which will be shown on screen through treatment in electronic circuits and microprocessor.

3.2 Measurement Precautions

- (1)Please check the device before using, and confirm that it can work normally.
- (2)The SpO₂ sensor and photoelectric receiving tube should be arranged in a way with the subject's arteriole in a position there between.

connected to a computer, the medical electrical systems should be in accordance with IEC60601-1-1.

(19)It is not recommended to use the equipment in electric knife and other high frequency interference environment.

(20)Children are not recommended to operate the equipment directly and they must operate the equipment under adult supervision.

(21)Keep the application environment free of dust, vibration, corrosion or combustible substances, and avoid temperature and humidity too high or too low.

(22)Turn off the equipment immediately if water splash the equipment or there is condensate in the equipment.

(23)Do not use the equipment immediately when it is transferred from cold environment to warm and humid climate.

(24)There is a visual indicator of low battery capacity in the equipment. Replace batteries in time when the low battery indicator appears.

(25)High temperature or high pressure steam disinfection of the Pulse Oximeter is not permitted.

(26)Do not have the Pulse Oximeter immersed in liquid. When it needs cleaning, please wipe itsurface with medical alcohol by soft material. Do not spray any liquid on the device directly.

(27)When cleaning the device with water, the temperature should be lower than 60°C.

(28)As to the fingers which are too thin or too cold, it would probably affect the normal measure of the

(3)The SpO₂ sensor should not be used at a location or limb tied with arterial canal or blood pressurecut or receiving intravenous injection.

(4)Do not fix the SpO₂ sensor with adhesive or else it may result in venous pulsation and inaccurate.

(5)Put the proper finger into the Pulse Oximeter. The measured value may be affected if the too cold or thin.

(6)The finger shall be inserted to the finger clip cavity deeply enough, and otherwise, the measured value may be inaccurate. Ensure the light between the photoelectric receiving tube and the luminescent tube of the Pulse Oximeter shall pass through the small artery bed of the finger of the subject.

(7)Keep the place where the light path passes free from light barriers such as adhesive tape to ensure the accuracy of parameter measurement.

(8)The product is suitable for children and adults with finger thickness of 7mm~25.4mm.

(9)The measurement will be affected by excessively strong ambient light, including but not limited to fluorescent lamp, double Ruby lamp, infrared heater, direct sunlight, etc.

(10)The measurement accuracy will be affected if the subject performs strenuous activities.

(11)The interference of electrical surgical equipment may affect the measurement accuracy.

(12)The nail on the test site shall be free from nail polish and other cosmetics. To ensure measurement accuracy, remove nail polish or other cosmetics if any.

Please take care when you insert the batteries for the improper insertion may damage the device.

6.3Mounting the Hanging Rope



Step 1. Put the end of the rope through the hole.

Step 2. Put another end of the rope through the first one and then tighten it.

7.OPERATING GUIDE

1.Insert the two batteries properly to the direction, and then replace the cover.

2.Open the clip as shown in Figure 8.

3.Let the testee finger put into the rubber cushions of the clip (make sure the finger is in the right position), and then clip the finger.

4.Press the button once on front panel.

5.Do not shake the finger and keep testee ase during the process.Meanwhile, human body is not recommended in movement status.


11.TECHNICAL SPECIFICATIONS

SpO ₂ Parameter Specification
Measuring range
Pulse Parameter Specification
Measuring range
Battery Requirement
Two batteries can work continually for 20 hours
Dimensions and Weight
Dimensions
Weight

Display Information	Display Mode
Pulse rate (PR)	Digital display
Pulse rate (PR)	Digital display
Pulse rate (Ba-graph)	Ba-graph display
Pulse wave	Waveform

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This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

RF warning for Portable device: The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.