User Manual

EVS 4343W / EVS 4343WP

EVS 3643W / EVS 3643WP

DRTECH

Digital Radiography Technologies

Manual Information

Copyright

Information contained in this manual is confidential and the property of DRTECH Corporation., which is only provided to the authorized customers of DRTECH Corporation. The purpose of the manual is to assist the users to comprehensively utilize the product of DRTECH Corporation. All information contained in the manual may not be altered or distributed without prior written consent of DRTECH Corporation., and may not be disclosed to unauthorized personnel. Information contained in the manual may be subject to change without prior notice.

Copyright © 2019 DRTECH Corporation. All rights reserved.

To Customers

Thank you for purchasing the DRTECH Radiography EVS4343W / EVS3643W / EVS 4343WP / EVS 3643WP(hereinafter this Product). This User's Manual explains how to use the detector, x-ray interface unit, and other peripheral equipment. Before using this product, be sure to read this manual thoroughly in order to utilize it more effectively. Also, read the Operation Manual for EVS Calibration and configuration Software (hereinafter ECali1).

Important information on usage and management of equipment

- Only a physician or legally certified operator should use this product.
- The equipment should be maintained in a safe and operable condition by maintenance personal.
- Use only computers and image display monitors complying with IEC 60601-1 or IEC 60950-1 and under a system configuration complying with IEC 60601-1. For details, consult your sales representative or local DRTECH dealer.
- Use only the dedicated cables. Do not use any cables other than those supplied with this product.

Disclaimer

- In no event shall DRTECH be liable for any damage or loss arising from fire, earthquake, any action or accident by a third party, any intentional negligent action by users, any trial usage, or other usage under abnormal conditions.
- Roentgenography, image processing, image reading, and image data storage must be performed in accordance with the laws of the country or region in which the product is being used. The user is responsible for maintaining privacy of image data.
- In no event shall DRTECH be liable for personal physical harm or property damage that is sustained, under any conditions when the instructions are not followed or the product is misused.
- It is the responsibility of attending physicians to provide medical care services. DRTECH will not be liable for fault diagnosis.

- In no event shall DRTECH be liable for direct or indirect consequential damages arising from the use or unavailability of this product. DRTECH shall not be liable for loss of image data for any reason.
- In no event shall DRTECH be liable for any damage arising from moving, alteration, inspection or repair by a person other than authorized service engineers.
- Specifications, compositions, and appearance of this product may change without prior notice.

Please contact us if there are any problems while using the product.

DRTECH Corporation

Customer Support Team

Website: <http://www.drtech.co.kr>

E-mail: drtech@drtech.co.kr

Tel: +82-31-779-7448

Fax: +82-31-779-7490

Definition of Symbols

Symbols	Description
	Non-ionized radiation
	The Waste Electrical and Electronic Equipment Regulations indicate separate collection for electrical and electronic equipment.
	Protective Earth (Ground)
	Direct Current
	Alternating Current
	Equipotentiality
	Power Off
	Power On
	Read and understand all instructions and warning labels in the product documentation before using the equipment. Keep this manual for future reference.
	Product contains specific materials that are suitable for recycling.
	Should be treated with care because if mistreated it might explode.
	Keep away from fire and flames.
	Heavy loading is prohibited.
	Authorized representative in the European Community
	Serial Number
	Country of manufacture (The date of manufacture may be added adjacent to this symbol.)
	Manufacturer
	For U.S.A standards Federal law restricts this device to sale by or on the order of a licensed practitioner
	For European Union (EEC Countries) This mark shows that product is in compliance with the essential requirements and other relevant provisions of Directive 93/42/EEC and "1639" shows the notified body number for MDD.
	Certification mark that indicates the product complies with ANSI/AAMI ES60601-1:2005/(R)2012 and A1:2012,, C1:2009/(R)2012 and A2:2010/(R)2012 and CAN/CSA C22.2 No.601.1, that specifies protection against fire, electric shock, and mechanical hazards.

Symbols	Description
	General Warning
	For Brazil It indicates that the product is in compliance with the essential requirements of Brazil National Telecommunications Agency.
	For Republic of Korea It indicates that the product is in compliance with Korea's product safety requirements for electrical and electronic equipment.
	Indicates the item is a medical device
	Indicates a carrier that contains unique device identifier information. (This symbol identifies the UDI carrier, including the AIDC and human readable information.)

TABLE OF CONTENTS

MANUAL INFORMATION	2
1. SAFETY AND REGULATORY INFORMATION	8
1.1 MEDICAL EQUIPMENT SAFETY STANDARDS	19
1.1.1 Medical Equipment Classification.....	19
1.1.2 Product Safety Standards.....	19
1.2 DECLARATION OF CONFORMITY	21
1.2.1 FCC Compliance Statement.....	21
1.2.2 IC Compliance Statement.....	22
1.3 LABELS AND MARKINGS ON THE EQUIPMENT	23
1.3.1 Detector.....	23
1.3.2 Other components	25
1.3.2.1 Battery Charger	25
1.3.2.2 Battery pack.....	25
1.3.2.3 EVS-WPCS	26
1.4 PRECAUTIONS REGARDING ELECTROMAGNETIC WAVES	26
2. DETECTOR SPECIFICATIONS	29
2.1 MAIN SPECIFICATIONS	29
2.1.1 Detector.....	29
2.1.2 Battery Charger.....	34
2.1.3 Battery Package	35
2.1.4 EVS-WPCS.....	36
2.2 CHARACTERISTICS	37
2.3 INTENDED USE	38
3. OPERATING PROCEDURE.....	40
3.1 NOTES FOR USING THE EQUIPMENT	40
3.2 SYSTEM CONFIGURATION.....	42
3.3 PARTS NAME AND FUNCTIONS.....	43
3.3.1 Product Components	43
3.3.1.1 Detector Components.....	44
3.3.1.2 EVS-CPDG Components.....	46
3.3.1.4 Battery Charger Components.....	47
3.3.1.5 Battery Pack Components.....	47
3.3.1.6 Charging Battery Pack.....	48
3.3.1.7 EVS-WPCS Components.....	51
3.3.2 Additional Function.....	51
3.3.2.1 AED Mode.....	52

3.3.2.2 Recommendation of setting AED Sensing Area.....	54
3.4 ENVIRONMENTAL REQUIREMENTS	55
3.4.1 Specifications of Workstation (Minimum and Recommended PC Requirements)	55
3.4.2 Recommended Specifications of Grid.....	56
3.4.3 Recommended Exposure Condition.....	56
3.5 BEFORE USING THE DETECTOR.....	57
3.6 OPERATING THE DETECTOR	58
3.6.1 Genral Workflow.....	58
3.6.2 Operating the detector	59
3.6.2.1 Powering on the detector.....	59
3.6.2.2 Register the detector and make connection to the EVS control system	59
3.6.2.3 Conduct examination	60
3.7 POWERING OFF THE DETECTOR	62
3.8 DETECTOR STATUS LIST	62
3.8.1 Detector Status List.....	62
3.8.2 Power Mode Indicator.....	63
3.8.3 Signal Strength Indicator.....	63
3.8.4 Battery Remains Indicator.....	64
3.9 ADDITIONAL FUNCTION.....	65
3.9.1 Lossless AED/AWC Mode.....	65
3.9.1.1 Recommendation of setting AED Sensing Area.....	67
3.10 PRODUCT INSTALLATION.....	68
3.10.1 Connecting Device.....	68
3.10.1.1 Operating AP	68
4. TROUBLE SHOOTING	70
5. MAINTENANCE	71
5.1 DAILY INSPECTION AND MAINTENANCE	71
5.2 REGULAR INSPECTION AND MAINTENANCE	71
5.3 IRREGULAR INSPECTION AND MAINTENANCE	71
6. PACKAGING	72
7. REVISION HISTORY	73

1. Safety and Regulatory Information

Intended Use

The EVS 4343W / EVS 3643W / EVS 4343WP / EVS 3643WP Digital X-ray detector is designed as a digital imaging solution, for providing radiographic diagnosis of the human anatomy. This device is intended to replace film or screen based radiographic systems in all general purpose diagnostic procedures. This device is not intended for mammography applications.

Safety Notices

Follow the safety instructions in this manual along with the following warning and cautions signs. Ignoring instructions, warnings, or cautions while handling the product may result in serious injury, accident, or product damage. To avoid any accident causing personal injury or product damage, be sure to read this manual carefully before using the product.

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Installation and Environment of Use

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- **Because the equipment cable is long, take care that cables do not become tangled during use. Also, be careful not to get your feet caught in the cable.**

Otherwise, it may cause a malfunction of the equipment or the injury of the user due to tripping over the cable.

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
- **Non-medical equipment such as the battery chargers, access point devices and wireless charging systems (EVS-WPCS) cannot be used in patient's vicinity.**

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WARNING

AVERTISSEMENT

- Do not use the battery pack as a power source for equipment other than EVS 4343W / EVS 3643W / EVS 4343WP / EVS 3643WP detector. Be sure to use the battery pack dedicated to the EVS 4343W / EVS 3643W / EVS 4343WP / EVS 3643WP detector.
- The battery charger is designed for the dedicated battery pack. Do not use the battery charger other than the dedicated one. Ignoring this may result in a battery explosion or a battery leak causing fire or electrical shock.
- Do not operate the battery charger using any type of power supply, other than the one indicated on the rating label.
- Do not handle the product with wet hands.
- Do not use the battery if battery connector is wet or moist.
- Do not attempt to disassemble, alter, or apply heat to the product.
- Avoid dropping or subjecting the product to severe impacts. To avoid the risk of injury, do not touch the internal parts of the battery if it has been cracked.
- Stop using the battery pack immediately if it emits smoke, a strange smell, or otherwise behaves abnormally.
- Do not let the battery pack and battery charger come into contact with water or other liquids and do not allow them to get wet.
- Do not clean with substances containing organic solvents such as alcohol, benzene, thinner, or other chemicals. Otherwise, fire or electrical shock may occur.
- Do not allow dirt or metal objects (such as hair pins, clips, staples or keys) to contact the terminals. Otherwise, battery explosion or leakage of electrolyte may occur, resulting in fire, injury or pollution of surrounding area. If the battery leaks and the electrolytes come into contact with your eyes, mouth, skin or clothing, immediately wash it away with running water and seek medical attention.
- Do not leave, store, or place the product in a location near heat sources, or in a place subject to direct sunlight, high temperature, high humidity, excessive dust, or mechanical shock. Otherwise, battery leakage, overheating or damage to the product may occur, resulting in electrical shock, burns, injury or fire.
- Do not attempt to use a battery pack that has deteriorated. Using a battery pack that has exceeded its life cycle may lead to overheating, fire or explosion.
- The Lithium ion/polymer battery is recyclable.
- Battery slowly discharges even if not in use.
- The useful battery life will be expired if it discharges immediately after being fully charged. You can purchase an optional battery pack to replace an exhausted one.
- The battery pack is a consumable item. If a fully charged battery is consumed quickly, use a new and fully charged battery pack.
- Be sure to charge the battery periodically (once a year) if it is not used for an extended period of time. The battery pack cannot be charged if it has been over discharged.
- Before discarding the battery pack, cover the terminals with adhesive tape

	or other insulators. Contact with other metal materials may cause fire or explosion. • Remove battery pack from equipment which is unused for a few minutes to prevent discharge.
--	--

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- **Have the patient take a fixed posture and do not let patient touch the parts unnecessarily.**

If a patient touches connectors or switches, it may result in electric shock or malfunction of the detector.

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- **Do not spill liquid or chemicals onto the equipment. Especially in cases that the patient is injured, which requires the equipment to come in contact with blood or body fluids, protect the equipment with a disposable cover as necessary.**

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- **Turn OFF the power to each piece of equipment for safety when not being used.**

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- **Handle the equipment carefully.**

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
- **Do not submerge the equipment in water.**

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- **The internal image sensor may be damaged if something hits against it, or if it is dropped, or receives a strong jolt.**

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- **Be sure to use the detector on a flat surface to ensure that the detector is not bended. Otherwise, the internal image sensor may be damaged. Be sure to securely hold the detector while using it in upright position. Otherwise, the detector may fall over, resulting in personal injury, or may flip over, resulting in damage to the inner components.**

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WARNING

AVERTISSEMENT

- Health care professionals and hospital administrators take special care in reducing radiation exposure to pediatric patients by following these steps:
- Discuss the rationale for the examination with the patient and/or parent to ensure a clear understanding of benefits and risks.
- Reduce the number of inappropriate referrals (i.e., justify X-ray imaging exams) by:
 - determining if the examination is needed to answer a clinical question,
 - considering alternate exams that use less or no radiation exposure, such as ultrasound or MRI, if appropriate, and
 - checking the patient's medical imaging history to avoid duplicate exams.
- Use the pediatric protocols or technique charts included as following table:

Table 1. Age categories

NAME	DEFINITION	FDA CODE
NEONATES	NEWBORNS UP TO ONE MONTH	NEO
INFANTS	ONE MONTH TO TWO YEARS	INF
CHILDREN	TWO YEARS TO TWELVE YEARS	CHI
ADOLESCENTS	TWELVE YEARS TO SIXTEEN YEARS	ADO
OTHER	OTHER AGE GROUP STUDIED	OTH

Table 2. Pediatric protocols

Age categories	NEONATES		INFANTS		CHILDREN		ADOLESCENTS		OTHER		Collimation size	Protection
Bodypart	kVp	mAs	kVp	mAs	kVp	mAs	kVp	mAs	kVp	mAs		
Skull AP	60	10	70	10	75	12	70	16	70	20	8 x 10 or 10 x 12	chest, abdomen
Skull LAT	60	10	65	10	70	10	70	16	70	20	8 x 10 or 10 x 12	chest, abdomen
Chest PA	70	1	90	2	90	2	90	4	90	4	10 x 12 or 14 x 14	abdomen, gonadal
Chest LAT	80	2	90	2.5	95	2.5	100	4	100	4	10 x 12 or 14 x 14	abdomen, gonadal
Abdomen	50	4	60	6.3	70	6.3	75	20	75	25	10 x 12 or 14 x 14	thyroid, gonadal
Pelvis	50	4	60	6.3	70	6.3	75	20	75	25	10 x 12 or 14 x 14	gonadal

Please refer to the FDA website for details on pediatric information.

[http://www.fda.gov/Radiation-](http://www.fda.gov/Radiation-EmittingProducts/RadiationEmittingProductsandProcedures/MedicalImaging/ucm298899.htm)

[EmittingProducts/RadiationEmittingProductsandProcedures/MedicalImaging/ucm298899.htm](http://www.fda.gov/Radiation-EmittingProducts/RadiationEmittingProductsandProcedures/MedicalImaging/ucm298899.htm)

1.1 Medical Equipment Safety Standards

1.1.1 Medical Equipment Classification

Type of protection against electrical shock	Class I ME Equipment
Degree of protection against electrical shock	Type B Applied Parts (Applied Part: Detector panel)
Degree of protection against ingress of water	IP54
Mode of operation	Continuous Operation
Flammable anesthetics	Not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide

1.1.2 Product Safety Standards

USA and Canada

ANSI/AAMI ES60601-1:2005/(R)2012 CAN/CSA-C22.2 No. 60601-1:14	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2 Ed.4: 2014	Medical electrical equipment-Part 1-2:Collateral standard: Electromagnetic compatibility-Requirements and tests

European Union

IEC 60601-1:2005 + C1:2006 + C2:2007 + A1:2012	Medical electrical equipment – Part1 : General requirements for Safety
EN 60601-1-2:2015	Medical electrical equipment–Part 1-2 : Collateral standard: Electromagnetic compatibility-Requirements and tests
IEC 60601-1-6 : 2010 + A1:2013	Medical electrical equipment –Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability
IEC 62304:2006 + A1:2015	Medical device software–Software life cycle processes

IEC 62366:2007 +A1:2014	Medical device – Application of usability engineering to medical devices
EN ISO 14971:2012	Medical device – Application of risk management to medical devices

The EVS 4343W / EVS 4343WP / EVS 3643W / EVS 3643WP may be operated in:

	BE	BG	CZ	DK	DE	EE	IE	EL
	ES	FR	HR	IT	CY	LV	LT	LU
	HU	MT	NL	AT	PL	PT	RO	SI
	SK	FI	SE	UK				

Abbreviations

Belgium (BE), Bulgaria (BG), Czech Republic (CZ), Denmark (DK), Germany (DE), Estonia(EE), Ireland(IE), Greece(EL), Spain(ES), France(FR), Croatia(HR), Italy(IT), Cyprus(CY), Latvia(LV), Lithuania(LT), Luxembourg(LU), Hungary(HU), Malta(MT), Netherland(NL), Austria(AT), Poland(PL), Portugal(PT), Romania(RO), Slovenia(SI), Slovakia(SK), Finland(FI), Sweden(SE) and United Kingdom(UK)

 WARNING	This device is restricted to indoor use in the 5150-5350MHz band.
--	---

1.2 Declaration of Conformity

1.2.1 FCC Compliance Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

FCC ID

FCC ID	Model	Remark
RNH-EVS4343W	EVS 4343W / EVS 4343WP EVS 3643W / EVS 3643WP	Component Modular Information Manufacturer: Realtek Model Name: RTL8812BU

WARNING

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

1.2.2 IC Compliance Statement

This Class B digital apparatus complies with Canadian ICES-003.

This device complies with Industry Canada license-exempt RSS standard(s).

Operation is subject to the following two conditions:

- (1) this device may not cause interference, and
- (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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For product available in the USA/Canada market, only channel 1~11 can be operated. Selection of other channels is not possible.

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This device and its antenna(s) must not be co-located or operation in conjunction with any other antenna or transmitter.

The device could automatically discontinue transmission in case of absence of information to transmit, or operational failure. Note that this is not intended to prohibit transmission of control or signaling information or the use of repetitive codes where required by the technology.

The device for the band 5150-5250 MHz is only for indoor usage to reduce potential for harmful interference to co-channel mobile satellite systems.

High-power radars are allocated as primary users (meaning they have priority) of the bands 5250-5350 MHz and 5650-5850 MHz and these radars could cause interference and/or damage to LE-LAN devices. 5600-5650MHz is restricted to be used at Canada.

1.3 Labels and Markings on the Equipment

The EVS 4343W / EVS 3643W / EVS 4343WP / EVS 3643WP detector and other components have labels and markings on them. Their contents and locations are indicated in the following chapters.

Any label removal or evidence of attempting to remove the label will void the warranty of the product.

1.3.1 Detector


Figure1. EVS 4343W & EVS 3643W Markings

This mark indicates that this is Type B Applied Part according to ANSI/AAMI ES60601-1 (2005) + AMD 1 (2012), CAN/CSA-C22.2 No. 60601-1 (2014) and EN 60601-1

This Mark indicates that this equipment must be handled with care.


Figure2. EVS 4343W / EVS 4343WP Labels



Figure3. EVS 3643W / EVS 3643WP Labels



Caution: Do not jolt or apply excessive load.

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
1.4 Precautions regarding electromagnetic waves

The Product is intended for use in the electromagnetic environment specified below. The user of the Product should assure that it is used in such an environment.

Emission Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR11	GROUP 1	The Product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electromagnetic equipment. The Product is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Compliance	

Electromagnetic immunity

The Product is intended for use in the electromagnetic environment specified below. The user of the Product should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic discharge(ESD) IEC 61000-4-2	$\pm(2, 4, 8)$ kV contact $\pm(2, 4, 8, 15)$ kV air	$\pm(2, 4, 8)$ kV contact $\pm(2, 4, 8, 15)$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines (100kHz PRF) ± 1 kV for input/output lines(100kHz PRF)	± 2 kV for power supply lines (100kHz PRF) ± 1 kV for input/output lines(100kHz PRF)	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to ground	± 1 kV line to line ± 2 kV line to ground	Mains power quality should be that of a typical commercial or hospital environment
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_T for 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T for 1 cycles 40 % U_T for 5 cycles 70% U_T for 25/30 cycles 0% U_T for 250/300 cycles	0% U_T for 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U_T for 1 cycles 40 % U_T for 5 cycles 70% U_T for 25/30 cycles 0% U_T for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EVS 4343W / EVS 4343WP / EVS 3643W / EVS 3643 WP requires continued operation during power mains interruptions, it is recommended that the EVS 4343W / EVS 4343WP / EVS 3643W / EVS 3643WP be powered from an uninterruptible power supply

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of atypical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3Vrms 150kHz-80MHz 6Vrms in ISM bands between 0.15 MHz and 80 MHz Amateur 80% AM at 1 kHz	3Vrms 150kHz-80MHz 6Vrms in ISM bands between 0.15 MHz and 80 MHz Amateur 80% AM at 1 kHz	The RF communications equipment should be used no closer to any part of the Product including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separations distance
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2700 MHz 80% AM AT 1 kHz	10 V/m 80 MHz to 2700 MHz 80% AM AT 1 kHz	$d=2.0\sqrt{P}$ 80 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation

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

- **Reliability (Lifetime Dose)**


More than 120Gy (206.9uGy/1 shot x200 shot/1day x 25days/Month x 12month/year x 10years = 124.14Gy)

2.1.2 Battery Charger


Item	Description
Model	EVS-BCS
Simultaneous Charging	Battery Pack 2 EA
Charging Time	3 hours
Rated Power Supply	DC +12 V, 6 A Max.
Dimensions (W x H x D)	180 mm x 255 mm x 35 mm
Weight	0.5 kg

2.1.3 Battery Package

Item	Description
Model	EVS-MBP-Y
Cell Type	Lithium Polymer
Number of Cells	2S1P (2series 1 Parallel)
Rated Power Supply	Output : DC +7.4 V
Lifetime	Approx. 500 cycles of use (complete charge/discharges 1 cycle)
Dimensions (W x H x D)	163 mm x 148 mm x 7 mm
Weight	0.24 kg

2.1.4 EVS-WPCS

Item		Description	Note
Model		EVS-WPCS	
Dimensions (W x H x D)		115.8 mm x 94.8 mm x 12 mm	
Weight		0.15 kg	
Charging Transceiver IC		Freescale MWTC1012	Medium Power
Charging Receiver IC		Freescale MWPR1516	Medium Power
WPC Qi Specification		WPC MP-A2 Standard.	
Rx, Tx Distance	Coil to Coil	Typ. 4mm(±1mm)	Max. 8mm
	Center to Center	Typ. 4mm	Max. 8mm
Input Voltage		DC 12V (±10%) / over 3A or Adaptor(XP Power / AHM85PS12 / DC 12V, 7.08A)	Tx Module Input
Output Power		DC 10V /1.4A (14W)	Rx Module Output
Standby Current/Power		Typ. 27.22mA / 326.6mW	
Max Power Efficiency		83%	
Ambient Temperature	Storage Temperature	-20℃ ~ 85℃	
	Operating Temperature	0℃ ~ 35℃	

2.2 Characteristics

- **Typical Patient Doses**

- Typical patient doses are equivalent to 500-1000 speed film/screen systems.

- **Sensitometric Characteristics and Dynamic Range**

- EVS 4343W/ EVS 4343WP / EVS 3643W / EVS 3643WP wireless responds linearly against the exposure range for 500-1000 speed film/screen where it can depict the clinical information. It means that products fully cover a dynamic range of 0.2-20 μGy at least.

- **User Friendly OLED Display**

- With its clear status indication, EVS 4343W / EVS 4343WP / EVS 3643W / EVS 3643WP's OLED Display provides various information at a glance. This function enables users to see the current status of the detector and the number of images stored in it to provide them better and more differentiated customer experience.

2.3 Intended Use

- **Intended Patient Population**

Considerations		Requirement Description
Age		Paediatric to Geriatric
Gender		No special requirements
Weight		>20 kg
Height		Not relevant
Nationality		No special requirements
Health		Not relevant
Patient state	Patient is user	Patient is not user
	Patient is not user	Not relevant, unless patient is agitated

- **Intended Part of the Body or Type of Tissue Applied to or Interacted with**

Measuring site: Overall body except to breast

Condition: Skin

- **Target diagnosis group**

- The group of patients is all general-purpose diagnostic procedures (excluding fluoroscopic, angiographic, and mammographic applications)

- **Target disease**

- ICD-11 Code: XY9R Diagnosis confirmed by imaging
- All general-purpose diagnostic procedures

- **Clinical benefits**

- non-invasive and painless diagnosis of disease and monitoring of therapy
- support of medical and surgical treatment planning
- Improvement of psychological
- Improvement of discomfort

- **Intended User Profile (Operator Profile)**

Considerations		Requirement Description
Education	Minimum	▪ A professional with specialist knowledge/ qualifications who has fully understood the product and the contents of this document. (Such as a doctor or radiological technologist)
Knowledge	Minimum	▪ Read and understand Western "Arabic numerals" written in Arial font ▪ Can distinguish human body ▪ Understands hygiene
Language understanding	Minimum	▪ Understand the manual in English

Considerations		Requirement Description
Experience	Minimum	Licensed physician or other qualified medical professional`
Permissible impairments	Minimum	- Mild reading vision impairment or corrected vision to log MAR 0.2 - Average degree of age-related short term memory impairment - Impaired by 40 % resulting in 60 % of normal hearing at 500 Hz to 2 kHz

 CAUTION MISE EN GARDE	U.S. Federal Law restricts this device to the sale by or on the order of a licensed physician or other qualified medical professional
--	---

- Intended Conditions of Use**

Considerations	Condition
Environment including hygienic requirements	<u>Operating conditions</u> - Temperature: +10 °C to +35 °C - Barometric Pressure: 700 hPa to 1060 hPa - Humidity: 20 % R.H. to 85 % R.H. <u>Storage and delivery conditions</u> - Temperature: -15 °C to +55 °C - Barometric Pressure: 500 hPa to 1060 hPa - Humidity: 10 % R.H. to 90 % R.H. Non-sterile Multiple patient use Less than ten minute contact Indoor use only Ambient luminance range: 100 lx to 1500 lx Viewing angle: normal to the display $\pm 20^\circ$
Frequency of use	1 day: 200 shot
Location	In hospital environment
Mobility	Portable ME equipment to be used on a patient

3. Operating Procedure

3.1 Notes for using the equipment

When using the equipment, take the following precautions. Otherwise, problems may occur and the equipment may not function properly.

System Diagnostic

The Ecali1 software generates a system diagnostic. Run Ecali1 software after installing the system and execute the program at least once a year. If an error occurs, report the detailed error information to DRTECH local dealer or distributor.

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The owner is responsible for ensuring that the system diagnostic is performed at least once a year. Do not try to use the system if an error is detected.

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Calibration

To ensure optimal performance of the system, it is important to verify that the system is calibrated.

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The owner is responsible for ensuring that the system calibration is performed after the system installation is completed or the system is repaired. Do not try to use the system if system calibration is not performed.

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Before exposure

- Make sure to check the equipment daily and confirm that it works properly.
- When moisture-laden air comes into contact with a cold surface, it will cause condensation to form on the equipment. In this case, wait until the condensation evaporates before performing an exposure.
- If the equipment is used while the condensation formed on it, problems may occur in the quality of the captured image.
- When an air-conditioner is used, be sure to raise/lower the temperature gradually so that temperature changes will not affect the product and prevent condensation.

During exposure

- Do not use the detector nearby devices generating a strong magnetic field. Doing so may produce image noise or artifacts.

Electric Shock Hazards

- To reduce the risk of electric shock, the system must be connected to an electrical ground.
- A three-phase AC power cable is supplied with this system to provide the proper electrical grounding. The power cable must be plugged into an UL-approved three-contact electrical outlet.
- Do not disassemble or modify the product as it may result in fire or an electric shock. There are no operator serviceable parts or adjustments inside the system. Only a trained and qualified person should be permitted access to the internal parts of the system.
- If an APPLIANCE COUPLER or Mains Plug or other separable plug is used, functional isolation device should be used for safety.

Disinfection and cleaning

- Wipe it with a dry cloth slightly dampened with a neutral detergent.
- Do not use solvents such as alcohol, thinner or benzene. Doing so may damage the surface of the equipment.
- The power should be switched off and equipment should be unplugged before cleaning.
- The exterior of the array can be cleaned with common hospital decontamination solutions including 5% of ButylCellosolve. A 0.55% Benzalkonium Chloride, 0.63% Sodium Hypochlorite, or 70% alcohol solution can also be used. To apply the cleaning solution, power down the system and disconnect from the power source, moisten cloth with the solution, and wipe the panel.
- The above procedure applies to all components.

Operating/storage environment

- Be sure to use and store the equipment under the conditions described below.

	Operation	Storage	Transit
Temperature	10 to 35 °C	-15 to 55 °C	-15 to 55 °C
Humidity	20 to 85% RH (Without Condensing)	10 to 90% RH (Without Condensing)	10 to 90% RH (Without Condensing)
Atmospheric pressure	700 to 1060 hPa	500 to 1060 hPa	500 to 1060 hPa

- Do not expose the equipment to high temperatures and/or high humidity. This may cause the product to malfunction.
- When not in use, keep the detector, handle unit, and grid in a designated location or in a location where they are safe and cannot fall down.

Notes on disposal

- Disposal of this product in an unlawful manner may have a negative impact on health and on the environment. When disposing of this product, therefore, be absolutely sure to follow the procedure which is in conformity with the laws and regulations applicable in your area.
- Refer to lifetime dose at "2.1 Reliability (Lifetime Dose)".

Handling the equipment

- The equipment must be handled with care to avoid personal injury or damage to the internal image sensor.

3.2 System Configuration

The Product is used in system configuration as illustrated below.


Figure 5. EVS W Series Single Detector Configuration

- **Do not connect other products to the remaining Lan port on the Ethernet Switching Hub. It may affect the image transmission of the detector.**
- **The following technical specifications represent recommended for the Ethernet Switching Hub.**
 - Supports 1000BASE-T or 1000BASE-TX
 - Supports full-duplex
 - Supports IEEE 802.3ab
 - More than 4 LAN Port Ethernet Hub

3.3 Parts Name and Functions

3.3.1 Product Components

Components	Model No.	Functions	Note
Flat Panel X-ray Detector	EVS 4343W EVS 4343WP EVS 3643W EVS 3643WP	Imaging device with Gigabit Ethernet signal interface to the operating PC	
Battery Charger (Incl. AC Power Adaptor, Power Cable)	EVS-BCS (Power Adaptor: AHM85PS12)	A device to charge a battery pack. Possible to charge two batteries packs simultaneously.	Power Cable: 2m
Battery Pack	EVS-MBP-Y	A device to supply power to the detector in Wireless Mode	2 EA
Wireless Charging System	EVS-WPCS	A device to wirelessly charge a battery pack that is attached to the detector	Optional Component
Power Adaptor (Incl. AC Power Cable)	MANGO60S-USB-PDA	A device to supply power to a Battery Charger	Optional Component Power Cable: 2m
Software CD	EConsole 1	Radiological Image Processing Software	Optional Component CE No. KR19/81826250

Table 3 EVS W Series Components Information

WARNING

The use of accessories and cables other than those specified, with the exception of **EVS W Series** accessories and cables sold by DRTECH Co., LTD. as replacement parts for internal components, may result in increased emissions or decreased immunity of the equipment. Accessory equipment connected to the analog and digital interfaces must be certified according to the respective IEC standards. All combinations of equipment must be in compliance with IEC 60601-1 system requirements. Any person who connects additional equipment to the signal input or signal output ports configures a medical system, and is therefore responsible for ensuring that the system complies with the requirements of the system standard IEC 60601-1. If in doubt, consult DRTECH technical support representative.

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3.3.1.1 Detector Components

The detector is designed to capture radiographic images. Captured images are transmitted to the image-capture computer using the wired data transfer.


Figure 6. Detector Components

A. Wireless Antenna: Transmit image data with wireless communication (IEEE802.11n).

- B. Battery Pack: Supply electrical power to the detector while using wireless communication
- C. Status Indicators
 - Power: Show power on/off status of the detector.
 - Ready: Show data communication status and ready status of the detector
 - Light ON: Detector is ready
 - Light OFF: Detector is busy
 - Link: Show detector's registration and connection status.
 - AP: Lamp indicating Wired/ Wireless mode (5 GHz)
- D. Display Button : Use to change OLED screen
- E. AP Button: Can select different wireless connection options for the detector.
(Connection options: Wireless using AP/ Wireless using detector's internal AP/ Portable mode)
- F. Power Button: Turn Detector on/off
- G. Connector: power supply through a USB Type-C cable or Power adaptor
- H. Debugging Cover
- I. WPCS Window: Wireless Charging Rx Window

 Simultaneous blinking of LED lamps indicate a system error


- Don't remove the item G and H. If the metal objects come into the case, it can lead to product malfunction. The operator should not make the patient contact inside of the access cover (G and H).

3.3.1.2 EVS-CPDG Components


Item	Description
Model	EVS-CPDG
USB PD Specification	Support PD Ver3.0
USB Version	USB 2.0
Ethernet	Support 100/1000BASE-T
Dimensions (W x H x D)	45.8 mm x 70.2 mm x 21.6 mm
Weight	0.1 kg

3.3.1.4 Battery Charger Components

- A. Power Indicator : Indicate the power on/off status..
- B. Charging Indicator : Indicate the charging status.
- C. Battery Compartment : Insert the battery pack to charge.
- D. DC Input : Connect the DC Adaptor to supply electrical power to the battery charger

3.3.1.5 Battery Pack Components

- A. Charging Indicator : Indicate the charging status
- B. Knob Latch : Rotate between on/off for battery swap


3.3.1.6 Charging Battery Pack

The battery pack supplies power to the detector. Be sure to use only the dedicated battery pack, and fully charge it before use.

- Connect the power adaptor to the DC Input port of the battery charger. The power LED light in blue indicates the direct current (DC) power is ON.
- Insert the battery pack into the battery charger. Charging starts automatically. The charge LED lights appear green when the battery pack is being charged. When the battery pack is completely charged, all levels of charge LEDs will illuminate.
- Gently pull the charged battery pack to remove it from the battery charger.
-


 WARNING AVERTISS-EMENT	Securely plug the power cord into the power source. If contact failure occurs, or if dust or metal objects come into contact with the exposed metal prongs of the plug, fire or electrical shock may occur.
 CAUTION MISE EN GARDE	Be sure to stop charging the battery pack when the charge LED lights appear in green longer than the specified charging time. Not doing so may result in battery pack overheating, smoke emission, battery explosion, or fire. You must use the power adaptor that is certified with IEC 60950 or IEC 60601-1.
 CAUTION MISE EN GARDE	Do not insert the battery pack into the charger in the opposite direction. If battery pack insertion failure occurs, the connector of battery pack is broken and the battery pack may not be charged or it may be damaged by electrical shock.
	Two batteries can be charged at the same time. It takes approximately two hours to fully charge a battery pack. The required charging time may vary depending on the temperature and remaining battery level.

How to Attach a Battery Pack

- 1) Align the arrows on the detector and battery pack.
- 2) Push down the battery pack.
- 3) Turn the battery lock knob 90 degrees clockwise.

How to Detach a Battery Pack


- 1) Turn the battery lock knob 90 degrees counter-clockwise.
- 2) Pull up the battery pack grabbing the knob.

 WARNING AVERTISS-EMENT	• Make sure to turn off the detector before detaching a battery pack. Press and hold the power button for about 2 seconds. All status LED lamps will be turned off when the detector is turned off. • Securely attach the battery into the detector or charger. If contact failure occurs, or if dust or metal objects come into contact with the exposed connector pins of the detector or charger, fire or electrical shock may occur.
--	---


How to Charge Battery Packs

Horizontal Direction (Attachment)


- Align the arrows on the charger and battery pack.
- Push down the battery pack.

Horizontal Direction (Detachment)


- Put the finger into the groove on the charger and grab the battery pack.
- Pull out the battery pack.

Vertical Direction (Attachment)

- Stand the battery pack upright to reveal the battery charged connector.
- Check the yellow-green arrows on the charger and bottom of the battery pack.
- Align the left and right side of battery pack to the charger.
- Push down the battery pack.

Vertical Direction (Detachment)

- Grab the battery pack.
- 2) Pull out the battery pack from the charger while holding the charger by the other hand.

3.3.1.7 EVS-WPCS Components


- ① EVS-WPCS Base
- ② RDA Window: IRDA communication window
- ③ Power Connection Cable
- ④ Debug Connector: For debugging
- ⑤ ID Switch: Device ID setting dip switch
- ⑥ Indicate LED Connector: For cradle (Optional)
- ⑦ Indicate LED Connector: For system (Optional)

3.3.2 Additional Function

Mode	Description
AED Mode	1. The detector detects actual amount of X-rays without any connection to the X-ray generator, and then performs image acquiring to the extent of image acquisition time and transits the image data. 2. No signal used (No need of connector interface cable) • You can use AED mode without connecting the generator with USB SW Box physically.

3.3.2.1 AED Mode

AED Mode is available for acquiring images without any connection to X-ray generator. Generator interface cable is not required.

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

3.3.2.2 Recommendation of setting AED Sensing Area

Standard Environment

We suggest the collimated area on detector is wider than 4cm X 8cm, and keep along the vertical direction as shown in below figures.

3.4 Environmental Requirements

3.4.1 Specifications of Workstation (Minimum and Recommended PC Requirements)

Item	Specification
Operating system	Windows 7 64 bit SP1 (Professional Edition or higher)
CPU	Intel Core i5 3470 or higher (or compatible CPU)
Memory	8GB or higher
Hard disk	1TB or higher
LAN card	Gigabit (Detector only) Intel® PRO 1000 Series (Gigabit LAN Card for network interface) Min. Requirements: 1Gbps, Jumbo Frames: 9K Receive Descriptors: 2K (higher than 1024) This is not dedicated to DICOM
Monitor	23 inch, 1920 x 1080 or higher
Optional disc drive	CD or DVD R/W

Table 4 Workstation

3.4.2 Recommended Specifications of Grid

Item	Recommended Specification
SID	100cm / 130cm / 150cm / 180 cm
Ratio	10 : 1
Frequency	103 line/inch, 120 line/inch, 215 line/inch
Inter spacer	Al

Table 5 Grid spec

3.4.3 Recommended Exposure Condition


Item	Recommended Condition
X-ray Energy Range	100cm / 130cm / 150cm / 180 cm
Reliability (Lifetime Dose)*	120Gy (206.9uGy/1 shot x200 shot/1day x 25days/Month x 12month/year x 10years = 124.14Gy)

Table 6 Exposure Condition

*Results obtained in a laboratory environment.

3.5 Before using the detector

 CAUTION MISE EN GARDE	Be sure to use only the specified power supply dedicated for the Product. Veillez à n'utiliser que l'alimentation appropriée du détecteur produit.
	Wireless Connection EVS W series transmits images and data by wireless communication. Battery pack should be installed in the detector to use it under the wireless configuration. Up to 2 battery packs can be charged simultaneously from a battery charger. Wired Connection Connect EVS W Series and PC with a tether interface and functional cable to make a wired configuration. Although the functional cable supplies power, a battery pack is needed to be installed in the detector. Data communications are faster than wireless connection. It is possible to charging a battery pack while using the detector.

3.6 Operating the detector

3.6.1 Genral Workflow

The following workflow indicates the procedures for acquiring a clinical image after startup of EConsole1 and other system equipments

- Preparing to use the detector

Attach a fully-charged battery pack to the detector.

¹⁾ : A procedure in order to register the detector to a specific digital radiography system

- Operating the detector

1. Turn on the detector

²⁾ : Network connection between the EVS 3643W wireless detector and the Ecali

2. Register¹⁾ the detector and make connection²⁾ to the EConsole1 power supply to the detector

³⁾ : • Selection of EVS 3643W wireless from the Exprimer series detector
• Selection of wireless/wired data transfer

3. Conduct Examination

• Select or register the patient information

Loop back procedure for each patient

• Select the protocol (selection³⁾ of the detector)

- Arrange the patient in the correct posture
- Position the X-ray generator to adjust the exposure field
- Check all the conditions

Loop back procedure for each body part

Check the captured images

- List the images
- Transmit the images

Sterilize the portion of the detector that has been in contact with a patient

• Conduct the next examination

- Ending use of the detector

1. Turn off the detector

2. Remove the battery pack

For details, refer to the operation manual or setup guide of the Ecali

Figure 8. General Workflow

3.6.2 Operating the detector

Lamp Type		Operation		Network
	Ready	Busy	AP	Link
Color	■Green	■Orange	■Blue	■Green

* The Network lamp is divided into two different lamps; AP (Blue) and LINK (Green) lamps.

3.6.2.1 Powering on the detector

- Before operating the detector, start up the EConsole1

Press and hold the POWER button
(approx. 1 second)

Operation lamp (Green) lights up

3.6.2.2 Register the detector and make connection to the EVS control system

Registration: AP(Blue) and LINK(Green) lamp blinks

- When the AP lamp is blinking once in 2 seconds, system is in wired mode.
- When the AP lamp is blinking twice in 2 seconds, system is in wireless mode (AP_1).
- When the AP lamp is blinking 3 times in 2 seconds, system is in wireless mode (AP_2).
- When the AP lamp is blinking 4 times in 2 seconds, system is in wireless mode (Detector AP)

**User can set value of AP_1, AP_2, in Ecali Program. Please refer to Ecali1 Operation Manual (Calibration tool).

Connection: Network connection between the internal wireless module of the detector and the wireless access point/EVS control system is secured automatically. The LINK (Green) lamp lights up when the detector is registered and the communication connection has been established.

- The LINK (Green) lamp does not light up when the detector is not registered or the communication connection is not established.
- When the READY (Green) and AP (Blue) lamp are blinking and LINK (Green) lamp does not light up, a communication error has occurred.

3.6.2.3 Conduct examination

For details about operation, refer to the User Manual for the EConsole1.

- ① **Select the patient information or protocols on the screen and start the examination.**

- The READY lamp color is changed from green to orange when the detector and EConsole1 change to exposure ready status.


- Arrange the patient in the correct posture and position the detector aligning it with the target body part.
- Position the X-ray generator to adjust the exposure field.
- Check all conditions before exposure.
- Make sure that Ready (Green) and LINK (Green) lamps are lit and AP lamp (Blue) is blinking. This means that the system is ready to start an examination

- A communication error has occurred when LINK lamp blinks.
- When the READY lamp (Orange) blinks slowly, the detector is in detector selection status (Sleep).
The detector enters detector selection status automatically when it has not been used for a certain period of time.

- ② **Press the exposure switch of the X-ray generator. Images captured with the detector are transmitted to the Ecali1and appear on the monitor.**

- Check the images on the monitor.
- If any uncompleted protocols remain, repeat the procedure(②).
- Choose the exposure mode before the shooting.

Mode	Description
AED	Auto Exposure Detection Mode

③ Click the button for ending the examination images are stored automatically.

- To conduct examination for another patient, repeat the step(③).

 IMPORTANT	• Sterilize the portion of the detector that has been in contact with a patient to prevent infection
---	--

3.7 Powering off the detector

- Turn off the detector

- ① Press the POWER button of the detector.
- ② All the LED lamps should be off.

- Select the POWER OFF button from Econsole1.
- Power plug off of POE adaptor from the outlet.

3.8 Detector status list

3.8.1 Detector Status List

Lamp Type		Operation		Network
	Ready	Busy	AP	Link
Color	■Green	■Orange	■ Blue	■Green
Power ON	O	X	X	X
During detector registration	X	O	☆	☆
Detector registration completed (1Sec.)	X	O	☆	O
During exposure preparation	X	O	☆	O
Ready status or performing an examination (Ready)	O	X	☆	O
During image data transmission	X	O	☆	O
Sleep Mode	X	O	☆	O
Deep Sleep Mode	X	O	☆	O
Power OFF	X	X	X	X

O : Light on

X : Lights off

☆ : Blinking

3.8.2 Power Mode Indicator

Display	Power Mode	Status	Required Actions
	Active		
	Sleep		Low power mode
	Deep sleep		Hibernation power mode only
	Power turned off or not linked	Disconnected communications	Power off

3.8.3 Signal Strength Indicator

Display	Signal Strength (communication stability)	Status	Required Actions
	Wireless, high (Stable)	Normal	
	Wireless, Normal (Stable)	Normal	
 	Wireless, Low (Unstable)	Unstable communication. Communication speed is lowered	Check whether there is any obstacle (e.g. your hands) between the wireless module and the wireless access point. If there is any obstacle, remove it. If the problem cannot be resolved, ask for consultation to your sales representative or local DRTECH dealer.
	No signal or No Link (Communication failed)	Disconnected communications	Confirm that detector and the access point are turned on. If the problem cannot be resolved, ask for consultation to your sales representative or local DRTECH dealer.
	Wired Link	Normal	External cable is connected.

- A signal strength indicator appears on the screen of the ECali1 computer. It shows the wireless communication level between the detector and ECali1.
- Keep the wireless communication level stable on capturing or transmitting images.

3.8.4 Battery Remains Indicator

Display	Status	Ext. Pwr	Required Actions
	Charge complement	Ext. cable & battery	
	Ext. cable charging	Ext. cable & battery	
	100%	Only battery	
	90~99%	Only battery	
	80~89%	Only battery	
	70~79%	Only battery	
	60~69%	Only battery	
	50~59%	Only battery	
	40~49%	Only battery	
	30~39%	Only battery	
	20~29%	Only battery	
	10~19%	Only battery	Warning message is popped up at the bottom-right. Recommend to change the battery.
	0~9%	Only battery	Warning message is popped up at the bottom-right. Change the battery before the battery is discharged.
	No Battery or Error	Unknown	Change the battery. If the problem cannot be resolved, ask for consultation to your sales representative or local DRTECH dealer.

3.9 Additional function

Mode	Description
AED/AWC Mode	3. The detector detects actual amount of X-rays without any connection to the X-ray generator. 4. No signal used (No need of connector interface cable)

3.9.1 Lossless AED/AWC Mode

AED Mode is available for acquiring images without any connection to X-ray generator. Generator interface cable is not required.

The acquired images are stored in the internal storage of the detector in order and the images can be transmitted to the PC after connecting the detector to the PC. The transmitted images in the detector are removed automatically.

Image processing and correction are available after connecting the detector to Viewer.

AWC(Auto Window Control) is a mode that automatically detects the end point of X-ray exposure and adjusts the EWT(Exposure Window Time) set before the examination.

If an X-ray is detected within a shorter time than the preset maximum Exposure Window Time, the Exposure Window is automatically closed and EWT is adjusted for faster operation of the detector.


Figure 9. AED Mode Configuration

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i In AWC mode, the MAX "Exposure Window Time" is set as a default in case the end of X-ray exposure is not detected.

3.9.1.1 Recommendation of setting AED Sensing Area

Standard Environment

We suggest the collimated area on detector is wider than 4cm X 8cm, and keep along the vertical direction as shown in figure 10.


Figure 10. Stand Environment with AED Mode

Table Environment

We suggest the collimated area on detector is wider than 4cm X 8cm, and keep along the horizontal direction as shown in figure 11.


Figure 11. Stand Environment with AED Mode

3.10 Product Installation

3.10.1 Connecting Device

This section describes how to connect the EVS 4343W / EVS 3643W system (Detector).

3.10.1.1 Operating AP

- 1) Connect the LAN cable to **ethernet** port(not internet) of AP, and the other to the LAN Card connector of workstation assigned for the data transfer.

- 2) Connect the power cable to the **power** port of AP to supply power.

- 3) Turn on **Wireless On/Off** switch and push the **Power** button.

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4. Trouble Shooting

When you encounter problems or error messages in the EVS Console while using this equipment, search the table below for the problem or error message and try the solutions.
If the problem persists, turn off the detector and consult your sales representative or local DRTECH dealer.
Please refer to the details of the problem or error messages.

Symptom	Cause/Error messages on the EVS Control Software	Remedy
Failed to turn on.	Not Installing battery pack well Dead battery pack Battery pack or Detector is broken	Install the battery pack. Charge the battery pack. Replace other battery packs and check the results. Replace other Detectors and check the results. Replace corresponding devices.
Ready LED(Orange) lamps of detector are blinking 3 time/sec	Internal hardware errors of the detector	Consult with service engineers of DRTECH
Consumption of a fully charged battery pack is fast.	Performance decrease caused by usage of long time. Usage of battery pack in low temperature environment.	Replace to new battery pack if the battery pack has been used for a long time. (Battery pack is a consumable.) Use battery pack in normal room temperature environment. Charging capacity of battery pack in a low temperature environment will decrease the capacity.
Battery pack or compartment for installation of battery pack is getting hot.	Battery pack failure. Detector failure.	Do not use the battery pack. Consult with service engineers of DRTECH.

5. Maintenance

In order to ensure that the equipment is used safely, be sure to inspect the equipment before use. If any problem is detected during the inspection and cannot be corrected, please contact your sales representative or local DRTECH dealer.

5.1 Daily Inspection and Maintenance

WARNING
AVERTISS-EMENT

For safety reasons, be sure to turn OFF the power to each piece of equipment before the following procedures. Otherwise, an electric shock may result.

- **Cable**
 - Ensure that cables are not damaged and cable jackets are not torn.
 - Ensure that the power cord plugs are securely connected to both the equipment AC inlet and the AC outlet.
- **Detector**
 - Ensure that there are no loose screws or broken parts.
 - Ensure that there is no dust or foreign substance on the external connector.
 - Ensure that there are no broken parts or short-circuits in the power supply connector.
- **After turning on the power**
 - Be sure to start the ECali1 before performing the inspection.
 - Perform exposure testing.


5.2 Regular Inspection and Maintenance

- **Monthly Inspection**
 - Ensure that there are no loose screws or broken parts.
 - Ensure that there is no dust or foreign substance on the external connector.
- **Yearly Inspection**
 - Perform a Performance Test or Self-diagnosis Test using a phantom or resolution chart, etc.

5.3 Irregular Inspection and Maintenance

- **Calibration**
 - Perform Calibration when exposure conditions have changed significantly. For details, refer to the Setup Guide for ECali1.

6. Packaging

7. Revision History

Revision	Date	Descriptions
00	2018.09.10	Initial Release
01	2019.09.09	Addition of EVS 3643WP / EVS 4343WPmodel in manual
02	2019.10.11	Update of regulatory information(FCC)
03	2019.10.14	Update of regulatory information(EU) and manual composition
04	2020.02.03	NB number change from 0120 to 1639
05	2020.02.13	Updating the RF output power values per frequency band
06	2020.08.11	Update of regulatory information(IC statement)
07	2021.04.09	- Address change of manufacturing site - Revision of the LABEL (Address, NRTL mark, symbol, warning text) - Revision of the RED Declaration of Conformity
08	2021.10.14	- Revision of the LABEL (KC mark, ANATEL mark) - Revision of "Definition of symbol" and typo
09	2022.02.08	Addition of the new LABEL for Medecom
10	2022.08.05	Update of FCC Compliance Statement
11	2022.12.10	- Revision of the LABEL (Address, NRTL mark, symbol, warning text) - Update of regulatory information(IC statement)

**DRTECH Corporation**

Suite No.1, 2 Floor / Suite No.2, 3 Floor, 29 Dunchon-daero541beon-gil,
Jungwon-gu, Seongnam-si, Gyeonggi-do, Republic of Korea

Website : <http://www.drtech.co.kr>

Customer Support Team

E-mail : drtech@drtech.co.kr

Tel : +82 31 7797448 / Fax : +82 31 7797490

**DRTECH Europe GmbH**

Am Kronberger Hang 2, 65824 Schwalbach am Taunus, Germany

D R T E C H