

深圳市晨北科技有限公司		料号 (Material No.)	XXX...
品牌 (Brand)	Etekcity	物料名称 (Material Name)	056
产品型号 (Model No.)	EFS-C671-WUST	材料+工艺 (Material+Process)	105g铜版纸, 6P, 风琴折, 单黑双面印刷
设计者 (Designer)	Coral	尺寸+公差 (Size+Tolerance)	148x210mm (LxW)
日期 (Date)	2025.01.03	版本号 (Version)	V0
备注 (Note)			
变更记录 (Revisions)			
		日期 (Date)	更改人 (By)
V0	初稿, 未受控	2025.01.03	Coral
A1	XXX...		
A2	XXX...		

READ AND SAVE THESE INSTRUCTIONS

Safety Information

Please read and follow all instructions and safety guidelines in this manual.

MEDICAL

WARNING: This scale uses a harmless electrical current to measure body fat. Consult a healthcare professional before using the smart scale if you use electronic medical equipment such as a pacemaker.

- This scale should not be used to diagnose or treat any medical condition.
- Scale data should not be used as a substitute for medical advice from a healthcare professional.

GENERAL USE

- Check the scale before using. **Do not** use the scale if it is damaged in any way.
- Do not** strike or bang on the scale.
- Do not** use the scale if there is a crack on the glass. While tempered glass is stronger and fractures more safely than ordinary glass, it may shatter unexpectedly after developing a slight crack. If you see a crack, contact **Customer Support**.
- The glass platform is slippery when wet. **Always** make sure the platform and your feet are dry before stepping on for measurement.

Not for commercial use.

Household use **only**.

Product working altitude: 0~3000M.

Product operating environment temperature: 5°C~50°C.

Product operating humidity: 35%~70%.

Pollution degree of the intended environment: 2 degree.

Ingress Protection Rating: IP20.

Leaving a battery in an extremely high temperature surrounding environment, such as sunshine, fire, or hot surface, may result in an explosion or the leakage of flammable liquid or gas;

A battery subjected to extremely low air pressure may result in an explosion or the leakage of flammable liquid or gas.

Do not mix new and old batteries and Do not mix alkaline, standard (carbon-zinc) or rechargeable (ni-cad, ni-mh, etc.) batteries.

BATTERY WARNING:

- Disposal of a battery into fire or a hot oven, or mechanically crushing in an explosion;

Federal Communication Commission Interference Statement – Part 15

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- This device may not cause harmful interference, and
- This device must accept any interference received, including interference that may cause undesired operation.

FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

FCC SDOC Supplier's Declaration Of Conformity

Etekcity Corporation hereby declares that this equipment is in compliance with the FCC Part 15 Subpart B. The declaration of conformity may be consulted in the support section of our website, accessible from www.etekcity.com

etekcity Smart Body Composition Scale

Quick Start Guide

1

To download the free VeSync app, scan the QR code or search "VeSync" in the Apple App Store® or Google Play Store.

2

Open the VeSync app. Log In or Sign Up.

3

Make sure your phone has Bluetooth® turned on.

4

Tap + in the VeSync app.

5

Step on the scale and select your smart scale from the device list.

6

Follow the in-app instructions to complete setup.

Note:

- Please make sure your firmware and software are up to date.
- The VeSync app is continually being improved and may change over time. If there are any differences, follow the in-app instructions.

Note:

- You can change the scale's in-app name and icon at any time in the smart fitness scale settings.
- When the scale display turns off, the scale's Bluetooth will disconnect and go into auto-off mode. However, when the scale display lights up and the app is opened to the device page, the scale will reconnect to Bluetooth.

Note:

- You can tap on the product categories to jump to different sections of the device list.
- Make sure to choose the right device model (you can find it on the back of the scale) or it will lead to a connection failure.

148 mm

210 mm

Display Messages

	Flashes when the scale needs to be charged.
	Flashes when Bluetooth is disconnected. Stays solid when connected.
	Weight exceeds 400 lb / 180 kg. Use a scale with a higher weight capacity.

Specifications

Name	Smart Fitness Scale
Model	EFS-C671-WUST
Weight Capacity	400 lb / 180 kg
Increments	0.1 lb / 0.05 kg
Weight Units	lb / kg
Platform	Tempered Glass
Battery	300mAh Li-ion Battery
Input	5.0V \equiv 500mA (USB Type-C Charging)
Dimensions	11.8 x 11.8 x 0.9 in / 30 x 30 x 2.4 cm
Operating Systems	iOS® 12.0+ or Android™ 6.0+

Using Your Smart Scale

- Place the scale on a hard, flat surface (not a carpet or mat).
- Step on the scale, then step off. Wait for the display to show **0.00** before weighing.
- Place dry, bare feet on the scale's conductive areas. Please stay on the scale until the screen shows your weight 3 times. If the scale is connected to the VeSync app, you can view additional body metrics in the app.

Note:

- If you are pregnant or have an internal medical device such as a pacemaker, turn on Zero-Current Mode in the VeSync app.
- To measure your baby or pet, tap in the VeSync app and follow the in-app instructions for weighing. Baby Mode will not measure other body metrics and will only measure weight and BMI.

Top

A Weighing Platform
B Display
C Conductive Areas

Bottom

D Reset Button
E USB-C Charging Port
F Unit Button

Tips

- If you're wearing shoes, the scale can only measure your weight. To measure other body metrics, such as body fat, step on the scale with **bare feet**.
- This scale uses a harmless electrical current to measure body fat. Consult a healthcare professional before using the smart scale if you use medical electronic equipment such as a pacemaker.
- If the scale has been moved or flipped upside down, it needs to be cleared.
 - Step on the scale until digits appear on the display, then step off.
 - The scale will show **0.0** when cleared.
- If your scale won't connect to the VeSync app:
 - Allow VeSync to access Bluetooth® in your phone settings.
 - Make sure to choose the right smart scale model (you can find it underneath the scale).
 - Make sure the scale is not connected to any other phone or app.
 - The scale should only be connected to your phone via the VeSync app. Please do not connect it in your phone's Bluetooth settings.
 - Turn off VPN if you are using one.
 - Download and open the newest version of the VeSync app.
 - Make sure your phone's operating system is iOS® 12.0+ or Android™ 6.0+.
 - For Android users, please allow the app to access nearby devices in your phone's settings.
- If you have other scale problems, try restarting the scale. Use a safety pin or paper clip to press the Restart button near the charging port to reset the scale. All saved data and configurations will be deleted.

Read the manual for further instructions, helpful tips, troubleshooting, and more.

Need More Help?

Scan here to view the manual for more information.

You can also type the following link in a web browser:
www.etekcity.com/pages/user-manuals

Quality Assurance

Your product comes with a 2-year limited warranty starting from the date of purchase.

- Scan the QR code or visit etekcity.com/warranty
- Enter your purchase info to receive exclusive news.
- Learn about our quality customer support and gain easy access to warranty information.

Customer Support

Etekcity Corporation

1775 Flight Way, Suite 150
Tustin, CA 92782
USA

Email: support@etekcity.com
Toll-Free: 1-855-686-3835

Support Hours

Mon-Fri
9:00 am-5:00 pm PST/PDT

**Please have your order invoice and order ID ready before contacting Customer Support.*

A0_254W_us

风琴折