


Care and cleaning

- Storage temperature range: It should be stored at room temperature between -25~+55°C, (-13°F~131°F) RH ≤85%, 700-1060hPa.
- Keep the unit dry and away from any liquids and direct sunlight.
- The probe should not be submerged into liquids.
- Please check the device if damaged once it falls. If you can't make sure of it, please send the complete device to the nearest retailer for re-calibration.
- Holding the thermometer too long may cause a higher ambient temperature reading of the probe. This could make the body temperature measurement lower than usual. Please try not to touch the thermometer probe part for avoiding the incorrect reading.

Installing Batteries

- Open the battery cover: use thumb to push battery cover out.
- Install or replace 2 "AAA" sized batteries in the battery compartment according to the indications inside the compartment.
- Replace the batteries in pairs. Remove batteries when the unit is not in use for extended periods of time.
- Place back the battery cover.



You need to replace the batteries when

- Low battery icon appears on the display.
- The ON/OFF/Memory button is pressed and nothing appears on the display.

Caution:

- Batteries are hazardous waste. Do not dispose them together with the household garbage.
- There are no user serviceable parts inside. Batteries or damage from old batteries are not covered by warranty.
- Use exclusively brand batteries. Always replace with new batteries together. Use batteries of the same brand and same type.

Troubleshooting

Error Message	Problem	Solution
EE	The system is not functioning properly.	Unload the battery, wait for 1 minute and re-power it. If the message reappears, contact the retailer for service.
E1	Exceeding operating temperature range	Allow the thermometer to rest in a room for at least 30 minutes at room temperature: for temple mode: 16~40°C (60.8~104°F). for object mode: 10~40°C (50~104°F).
H1	(1) In Temple mode: Temperature taken is higher than +42.2°C (108°F) (2) In Object mode: Temperature taken is higher than 99.9°C (199.9°F)	Please select the target within specifications. If a malfunction still exists, please contact the retailer for service.
Lo	(1) In Temple mode: Temperature taken is lower than +34°C (93.2°F) (2) In Object mode: Temperature taken is lower than 0°C (32°F)	Please select the target within specifications. If a malfunction still exists, please contact the retailer for service.

	Device can not be powered on to the ready stage.	Change with a new battery. If the message reappears, contact the retailer for service.
	Low battery: Battery icon is flashing, still can be measurable.	Suggest to replace the battery.
	Dead battery: Battery icon always enabled, no more measurements are possible.	Replace the battery

Specification

Temperature measurement range	Temple mode: 34~42.2°C (93.2~108°F) Object mode: 0~99.9°C (32~199.9°F)
Accuracy	For temple mode: ± 0.3°C (0.5°F) during 34~42.2°C (93.2~108°F). For object mode: ± 2°C or ± 3.6°F when measuring the object temperature between 0~99.9°C (or 32~199.9°F)
Operating condition	Temple mode: 16~40°C (60.8~104°F), RH ≤85%; 700~1060 hPa Object mode: 10~40°C (50~104°F), RH ≤85%; 700~1060 hPa
Storage and transportation condition	Temperature: -25~55°C (-13~131°F), RH ≤85%; 700~1060 hPa
Memory	9 set
Display resolution	0.1
Battery	AAA X 2
Weight (without battery)	160g (without battery)
Size	10.4cm(L) x 5.2cm(w) x 15.7cm(H)
Auto shutdown	60 sec.
Safety classification	⚠ Type BF equipment
IP Classification	IP22 Protection against water and dust

* Dispose of device and batteries according to local regulations.

Caution

- Only use the device once you have read and understood these instructions for use. Keep the instruction manual for use.
- The thermometer and patient needs to stay in the room in which the measurement is taken for at least 30 minutes before use.
- It is recommend that you always take the temperature in the same location, since temperature readings may vary according to the locations.
- It is recommended that three temperatures are taken and the highest one taken as the reading if:
 - New born infants
 - Children with a compromised immune system and for whom the presence or absence of fever is critical.
 - When the user is learning how to use the thermometer for the first time until he/she has familiarized himself/herself with the instrument and obtains consistent readings.
 - The reading is extremely low.
- Do not take a measurement while drinking, eating, sweating or exercising.
- Do not take measurement while or immediately after nursing.
- Children is not allowed to use the device. Medical products are not toys.
- Influences on body temperature
 - A person's individual metabolism
 - Age: Greater temperature fluctuations occur faster and more often in children. Normal body temperature decreases with age.
 - Clothing
 - Outside temperature
 - Time of day: Body temperature is lower in the morning and increases throughout the day towards evening.
 - Activities: Physical and, to a lesser extent, mental activities increase body temperature.
- This device is not shock-protected. Do not drop or expose to heavy shock this device.

- Do not bend the device.
- Do not disassemble or make modifications the device.
- Please do not dispose of the product in the household waste at the end of its useful life. Disposal can take place at your local retailer or at appropriate collection points provided in your country.
- Do not boil the probe.
- Do not use the device if it operates irregularly or any error message showed.
- Do not use thinner or benzene to clean the device.
- Wipe the device clean before storing.
- When take the device from storage at below or above 16~40°C(60.8~104°F), place it in 16~40°C(60.8~104°F) temperature range for at least 30 minutes before use.
- Remove the battery if the device will not be used for a long time.
- If this device is used according to the operation instruction, periodic re-calibration is not required. If you still have questions, please send the complete device to dealers.
- Do not use this device for diagnosis or treatment of any health problem on disease. Measurement results are for reference only. Contact your physician if you have or suspect any medical problems. Do not change your medications without the advice of your physician or healthcare professional.
- This device may not meet its performance specification if stored or used outside temperature and humidity ranges specified in specifications.
- Battery should not be charged or placed into extreme heat as it may explode.
- Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are inaccessible to children. If a battery has been swallowed, call a doctor immediately.
- Be sure to consult a doctor if you feel that your health is in poor condition.
- Do not judge your health only on the presence or absence of a fever.
- Do not service or maintain device while in use.
- The equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment, under 47 CFR 2.1093 paragraph (d)(2).

WARNING: The symbol on this product means that it's an electronic product and following the European directive 2012/19/EU the electronic products have to be dispose on your local recycling centre for safe treatment.

- Apple e il logo Apple sono marchi di Apple Inc., registrati negli Stati Uniti e in altri Paesi. App Store è un marchio di servizio di Apple Inc.
- Google Play e il logo di Google Play sono marchi di Google Inc.

Electromagnetic Compatibility Information

- This device needs to be installed and put into service in accordance with the information provided in the user manual.
- WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the HC700BT, including cables specified by the manufacturer. Otherwise, degradation of the performance of this device could result.

If higher IMMUNITY TEST LEVELS than those specified in Table 9 are used, the minimum separation distance may be lowered. Lower minimum separation distances shall be calculated using the equation specified in 8.10.

Manufacturer's declaration-electromagnetic immunity			
The HC700BT is intended for use in the electromagnetic environment specified below. The customer or the user of the HC700BT should assure that is used in such and environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms: 0,15 MHz – 80 MHz 6 Vrms: in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	Not applicable Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the HC700BT including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1,2 \sqrt{P}$, $d = 1,2 \sqrt{P}$ 80MHz to 800 MHz, $d = 2,3 \sqrt{P}$ 800MHz to 2,7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	

NOTE1: At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

More information on EMC compliance of the device may be obtained from GMC using the contacts shown in this manual.

Warranty Card

This instrument is covered by a 2 year guarantee from the date of purchase, batteries and accessories are not included. The guarantee is valid only on presentation of the guarantee card completed by the dealer confirming date of purchase or the receipt. Opening or altering the instrument invalidates the guarantee. The guarantee does not cover damage, accidents or non-compliance with the instruction manual. Please contact your local seller/dealer or GMC Inc.

Customer Name: _____

Address: _____

Telephone: _____

E-mail address: _____

Product Information

Date of purchase: _____

Store where purchased: _____

Issuance date: January 19, 2021

GMC Inc.
No. 686, Su Chu Rd.,
Chuzhou, Anhui, China



GMC_IB_HC700BT-
EN_AH_ver2104-
INDHC700BT000000XX