

Marking the CU RING with FCC ID:OSEALRU

Photographs of the exterior appearance of the CU Ring and its packaging are provided. The FCC identifier and a label will not be placed on the device itself for two reasons. First, the CU Ring is so small as to make placing the label on it is impracticable (Photograph 1). Second, the CU Ring is a disposable unit, will be always within its packaging including the time of purchase, until the time of use in a sterile environment. Having the label on the CU Ring would be of little use.

Photograph 1: AutoLap CU Ring general view




Photograph 2: AutoLap CU Ring back cover with a label



Therefore, consistent with FCC 47CFR § 2.925(f) the Applicant proposes as an alternative method of affixing the label, to affix it to the CU Ring packaging as shown in the photograph above and in the accompanying **User Manual**, provided as a separate exhibit to this application.

Sincerely,



Alon Shacham, VP R&D
Medical Surgery Technologies Ltd.