



Rayence Co.,Ltd.

14, Samsung 1-ro 1-gil, Hwaseong-si, Gyeonggi-do, Korea
Tel. +82-31-8015-6459 / Fax +82-31-8015-6394

Date: November 8, 2016

To: SIEMIC, INC.
775 Montague Expressway,
Milpitas, CA 95035
USA

Dear Sir/Madam,

Re: Letter of Similarity

We, Rayence Co.,Ltd., declare the differences for between original model and variant models.

-Product: Wireless digital X-ray flat panel detector
-Model name: 1012WCC
-Variant Model name: 1012WGC

There's different Sintillator type.

1012 WCC (Basic M/N) & 1012 WGC
CsI(Cesium Iodide/1012WCC) or Gd2O2S:Tb(Gadolinium Oxysulfide/1012WGC)

Sincerely Yours,

A handwritten signature in black ink, appearing to be "J8" or similar, written in a cursive style.

NAME : Namyoung, Ha
TITLE : Assistant Manager
E-MAIL: namyoung.ha@rayence.com
FAX: + 82-31-8015-6394
TEL: + 82-31-8015-6459
Rayence Co., Ltd.