




10 pcs

BioSensors
Ear tag BAT1





BioSensors Ear tag BAT1



10 pcs


BioSensors Ear tag BAT1




1. Read the RFID transponder



2. Check cow number




3. Add BA tag and read with NFC



4. Add backplate and tag to applicator



5. Dip applicator with backplate in disinfection



6. Tag the cow

Battery type: ER14250
Voltage: 3.6V
Frequency: 6240-6739.2Mhz
Art. no. 2150020580

This device complies with Part 15 of the FCC Rules.
Operation is subjected to the following conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.



FCC ID: UCSBAT1
IC: 6576A-BAT1

DeLaval International
Gustav De Laval's väg 15
SE-14751 TUMBA
Sweden



BioSensors
Ear tag BAT1

10 pcs



DeLaval

BioSensors
Ear tag BAT1

10 pcs






100 pcs

BioSensors
Ear tag
BAT1



BioSensors
Ear tag
BAT1



BioSensors
Ear tag
BAT1



100 pcs



BioSensors
Ear tag
BAT1

100 pcs



BioSensors
Ear tag
BAT1



BioSensors
Ear tag
BAT1

100 pcs

Battery type: ER14250
Voltage: 3.6V
Frequency: 6240-6739,2MHz
Art. no. 2150020581

This device complies with Part 15 of the FCC Rules.
Operation is subjected to the following conditions:
• This device may not cause harmful interference.
• This device must accept any interference received, including interference that may cause undesired operation.



FCC ID: UCSBAT1
IC: 6576A-BAT1
DeLaval International
Gustav De Laval's väg 15
SE-14751 TUMBA
Sweden



100 pcs

