

RF Exposure Exhibit

EUT Name: NeuroPace® Wand

Model No.: W-02

CFR 47 Part 2.1093

Prepared for:

Barbara Gibb
NeuroPace Inc.
1375 Shorebird Way
Mountain View, CA 94043
Tel: (650) 237-2700
Fac: (650) 237-2701

Prepared by:

TUV Rheinland of North America, Inc.
1279 Quarry Lane
Pleasanton, CA 94566
Tel: (925) 249-9123
Fax: (925) 249-9124
<http://www.tuv.com/>

Report/Issue Date: February 2, 2011
Report Number: 31051314.002

<i>RF Exposure Exhibit</i>	<i>1</i>
<i>1 RF Exposure</i>	<i>3</i>

1 RF Exposure

1.1 Exposure Requirements

The telemetry wand Model W-02 is used for programming and interrogating the implant device. It has no charging capability; neither the implant device can take the charge. According to CFR47 Part 2.1093(c) the Wand Model W-02 is categorically excluded from routine environmental evaluation for RF exposure. The exclusion is confirmed by the commission through TCB Permit But Ask Procedure. The referenced TCB PBA is attached as operational description exhibit.

For the current feeding the wand coil, please refer to page 6 of the W-02 Operational Description.

1.2 Test Procedure

Test is not required.

1.3 Evaluation

Na

1.4 Conclusion

SAR data is not required for FCC since the Wand Model W-02 is categorically excluded.