



Electronic Blood Pressure Monitor  
Model: HWA21  
Power: DC 3V(AA x2) 1.5W IP22  
Measuring Range: 0~295 mmHg  
FCC ID: 2ADYL-HWA21

This device complies with Part 15 of the FCC Rules.  
Operation is subject to the following two conditions:  
(1) this device may not cause harmful interference, and  
(2) this device must accept any interference received,  
including interference that may cause undesired operation.

Shenzhen Jumper Medical Equipment Co., Ltd.  
D Building, No. 71, Xintian Road, Fuyong  
Street, Baoan, Shenzhen, Guangdong, China

MedPath GmbH  
Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany

MD V1.0  
CE 0482  
2022-05  
Warning icons: crossed-out trash can, exclamation mark in triangle, person with heart, person with pulse line, and person with pulse line and heart.