


Appendix D

Photographs

1. SAR measurement System
2. Photographs of Tissue Simulate Liquid
3. Photographs of EUT test position
4. EUT Constructional Details

1. SAR measurement System



2. Photographs of Tissue Simulate Liquid

Photo 1: Tissue Simulant Liquid for HBBL600-10000MHz	NA
 A photograph showing a cross-section of a container filled with a yellowish, translucent liquid. A metal ruler is placed vertically next to the container, with the numbers 6, 7, 8, 9, and 10 visible. The liquid level is approximately at the 8.5 mark. The container has a dark blue or black outer shell and a lighter inner lining.	NA

3. Photographs of EUT test position

Photo 2: Left Cheek




Photo 3: Left Tilted



Photo 4: Right Cheek




Photo 5: Right Tilted



Photo 6: Front Side 10mm



Photo 7: Back Side 10mm



Photo 8: Left Side 10mm



Photo 9: Right Side 10mm



Photo 10: Top Side 10mm



Photo 11: Front Side 0mm



Photo 12: Back Side 0mm



Photo 13: Left Side 0mm



Photo 14: Right Side 0mm



Photo 15: Top Side 0mm



Photo 16: Front Side 14mm



Photo 17: Back Side 19mm



Photo 18: Top Side 14mm



4. EUT Constructional Details

