



Appendix A. FCC ID Label and Location

Product Label Sample with FCC Label information

Following is a sample copy of the label that will be placed on the rear cabinet of the product.
The FCC ID Label and compliance statement are marked in the product label.

| | |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------|
| Jawon medical CO., LTD. | |
| Automatic Blood Pressure Monitor | |
| FCC ID | : RIXFT-500 |
| Model NO | : FT-500L |
| POWER | : 230V~/50Hz(EC), AC 120V/60Hz(US) 50VA |
| Serial No | : N/A |
| <p>This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:</p> <ul style="list-style-type: none">(1) This device may not cause harmful interference, and(2) This device must accept any interference received, including interference that may cause undesired operation. | |
| <p>Manufacture : Jawon medical CO., LTD.</p> <p>Made in Korea</p> | |

Label Location

ETL Inc.
#584 Sangwhal-ri, Kanam-myon,
Yojiu-kun, Kyonggi-do,
469-880, Korea

Jawon medical CO., LTD.
Full Automated Blood Pressure Monitor
Model: FT-500L

Appendix A. FCC ID Label and Location

The label shown shall be permanently affixed at a conspicuous location on the device and be readily visible to the user at the time purchase.(Labeling requirements per 2.925)

